

Follow Up Renal Ablation

Siemens go.All

Application Examples: f/u a ablation	
Scan Protocol (See for Technical Factors)	Renal Mass
Oral Contrast	No
IV Contrast / Volume	Omnipaque 350 / P3T

This protocol is used to evaluate post ablation patients.

Patient Position: Patient lying in supine position, arms positioned comfortably above head and lower legs supported.

Scan Instructions & Range:

1. **Non-contrast & arterial phase** – scan kidneys only.
2. **90 second phase** – scan upper abdomen including entire kidneys.
3. **Delay images – 5-minute delay**

Recons & Reformations: Adjust FoV to fit body contour.